

MS Thesis (Graduation: 2019)

An Active Orthosis with Tenodesis Effect for Grasp Force Augmentation of Daily Life

Daegu Gyeongbuk Institute of Science and Technology (DGIST), Korea

Committee members:

1. Professor Jonghyun Kim
2. Dr. Taehun Kang
3. Professor Dongwon Yun

Background

✓ What is Tenodesis Effect?


- ✓ Wrist motion driven finger motion mechanism
 - Voluntary wrist extension
 - Minimal grasp force
- ✓ Wrist driven mechanical device assisting in making grasp for activities of daily living
 - C6-C7 spinal cord injury patients

✓ State of the art




- 3 motors
- Complex tendon routing
- Fixed grasp force (20N)
- Difficult to sanitize

➤ Robots with multiple motors/sensors



- 4 motors + IR + flex sensor
- Complex tendon routing
- Fixed grasp force (16N)
- Difficult to sanitize

➤ No way to adjust or tune grasp force





- Minimum grasp force
- Lower tasks accomplishments
- No force tuning

➤ Difficulty in sanitization






Experimental Setup

✓ Active 4-bar Mechanism



$$F = \frac{L_2 \cos(90 - \theta_3)}{L_2 \cos(90 - \theta_3) + L_1 \sin(\theta_1 + \theta_4)} \times \frac{1}{d} (T_m + T_w)$$




T_m = Constant motor torque, T_w = wrist torque



- Designed in SolidWorks
- Uses single MAXON BLDC motor with EPOS controller
- Mini gearbox with mini clutch mechanism for 90-degree force transmission
- Optimized linkage lengths through kinematics derivation
- Safety mechanism

Advantages

✓ Solved problems



- Generates up to 15N grasp force
- Accomplishes more than 75% of daily life activities
 - 10.4 N is required for 75% tasks
- Uses only 1 motor and no sensor
- Can operate in both active and passive mode
- Grasp force can be adjusted/tuned according to user's need
- Safety mechanism included to prevent excessive motor torque
- Easy to wear and sanitize

Other robots	My robot
Multiple motors & sensors	1 motor, 0 sensor
No force tuning	Force tuning
Sanitization issues	Easy sanitization